

## TrusMQ 16/48/96 Real-Time PCR System

### Product features:

- Tube/plate-flexible qPCR series that supports multi-color applications and multiplex protocols
- Rapid heating up and cooling down heating function
- Stand-alone operation with 7 or 10.4-inch touchscreen control
- Network operation with PC software control via direct connection of LAN
- Automatic recovery of the experiment when the power is on again after cutting off



| Category              | TrusMQ16 | TrusMQ48 | TrusMQ96 |
|-----------------------|----------|----------|----------|
| Throughput            | 16       | 48       | 96       |
| Fluorescence Channels | 2        | 4        | 6        |

### Reagent Kits :

| Cat. No                         | Product Name                                                                | Spec.    |
|---------------------------------|-----------------------------------------------------------------------------|----------|
| Lyophilization Detection Assays |                                                                             |          |
| RE002013                        | SARS-CoV-2 and Influenza A/B Virus Nucleic Acid Detection Kit               | 32 tests |
| RE002014                        | SARS-CoV-2 Nucleic Acid Detection Kit                                       | 32 tests |
| RE002015                        | Monkeypox Virus Nucleic Acid Real-Time PCR Detection Kit                    | 32 tests |
| RE002016                        | HPV Nucleic Acid Real-time PCR Detection Kit                                | 32 tests |
| RE002017                        | TALK/ROS1/MET Gene Fusions and Mutations Real-Time PCR Detection Kit        | 24 tests |
| RE002018                        | Mycobacterium Tuberculosis Complex Nucleic Acid Real-Time PCR Detection Kit | 32 tests |
| RE002019                        | Genitourinary pathogen nucleic acid real-time PCR detection kit             | 32 tests |
| RE002020                        | EGFR gene mutations real-time PCR detection kit                             | 24 tests |

| Cat. No                        | Product Name                                                                      | Spec.    |
|--------------------------------|-----------------------------------------------------------------------------------|----------|
| Liquid Format Detection Assays |                                                                                   |          |
| RE002005                       | DIAGNOVITAL <sup>®</sup> SARS-CoV-2 Real-Time PCR Kit                             | 50 tests |
| RE002006                       | DIAGNOVITAL <sup>®</sup> HS SARS-CoV-2 Real-Time PCR Kit v2.0                     | 50 tests |
| RE002007                       | DIAGNOVITAL <sup>®</sup> SARS-CoV-2 Multiplex Real-Time PCR Kit                   | 50 tests |
| RE002008                       | DIAGNOVITAL <sup>®</sup> INFLUENZA A/B SARS-CoV-2 Multiplex Real-Time PCR Kit     | 50 tests |
| RE002009                       | DIAGNOVITAL <sup>®</sup> INFLUENZA A/B SARS-CoV-2 RSV Multiplex Real-Time PCR Kit | 50 tests |
| RE002010                       | RTA <sup>®</sup> HBV DNA Real-Time PCR Kit                                        | 25 tests |
| RE002011                       | RTA <sup>®</sup> HCV Real-Time PCR Kit                                            | 25 tests |
| RE002012                       | RTA <sup>®</sup> HIV-I RNA Real-Time PCR Kit                                      | 25 tests |